

NEUTRON

ULTRASONICS

NEUTRON® technology, a guarantee of efficacy and safety. The perfectly-controlled vibrations and high-quality tips bring to treatments more preservation, efficacy and comfort.

Plus the widest range of tips on the market (80 different tips) with exclusive designs, alloys and coatings for clinical versatility. Exclusive tips in Pure Titanium for implant cleaning.



STEEL TIP
QUALITY:
THE STANDARD

NEUTRON Tips hardness, designed with a specific steel, is closest to enamel.

With our tips and linear displacement, teeth are not only cared for, but preserved

BUY 1 x P5XS B.LED or P5 B.LED or BOOSTER GET 1 x Handpiece FREE

PIEZOELECTRIC SURGERY IS POWERFUL ENOUGH
TO CARRY OUT EXTRACTIONS AT THE SAME SPEED¹
AS ROTARY INSTRUMENTS!

- Socket preservation⁹
- Safe for soft tissues
- Recommended for immediate implant placement
- No force needed
- 50% less pain & swelling²
- 98% less analgesics intake¹

1, 2, 9 References available on file



PIEZOTOME
CUBE

**BUY 1 x PIEZOTOME® CUBE (F50102),
GET 1 x Extraction II Tip Kit (F87546) or
1 x Universal Tip Kit (F87539) FREE**



SERVOTOME

High frequency electrosurge for orthodontics, periodontics and other oral procedures.

Efficient & Secure Surgery

Servotome totally respects biological tissues: the fine and even cut with no side effects allows good-quality healing.

30 Watts is Sufficient


Due to the technology evolution associated with the conductive bracelet, 30 watts is sufficient to perform all dental incision and coagulation procedures providing greater safety.



Experience the benefits of electrosurgery


Electrosection

Vaporizing the tissue with the active electrode makes it possible to cut quickly without damaging the edges of the incision.





Fulguration

Makes it possible to surface coagulate gingival tissue while protecting the underlying layer.



Coagulation

Concentrating energy at the surface of a massive electrode makes it possible to diffuse the heat to the surrounding tissue and produce instantaneous hemostasis.



BUY 1 x Servotome II, GET 1 x Electrode FREE

Acteon Ultrasonic Cleaners

HYGIENIC Aluminum control panel

SILENT PEHD cover

Stainless steel Structure / 2 years guarantee

Waterproof switch

ULTRASONICS CLEANER 5L - 280x180x260mm
- Usable volume 3.2L (PH500N)

ULTRASONICS CLEANER 5L WITH HEATER
20-60°C, 280x180x260mm
- Usable volume 3.2L with Heater (PH500CN)



Contact your preferred Acteon dealer for Special Ultrasonic Cleaner Pricing

IMPLANT CENTER 2

The alliance of technologies for safe and atraumatic bone surgeries

The alliance of the most advanced ultrasonic and rotary technologies ensures total independence in increasingly diverse clinical areas like: periodontics, endodontics, prophylaxis, restorative, pre-implant, implant surgeries.

Innovative, broad and robust surgical and non- surgical tips range



2 YEARS WARRANTY

BONUS

**BUY 1 x Implant Center 2® (F27203),
GET 1 x Cube LED Handpiece & Cord (F12816) FREE**

MINILED — SUPERCHARGED

Wavelength: 420 - 480 nm

Power: 2,000mW/cm²

Li-Ion battery - Handpieces:
anodized aluminium

Delivered with a Mini LED
handpiece, a battery,
a battery charger, a
power supply cord
with european plug,
a Ø7.5mm opalescent
guide and an oval
light shield



MINILED — ORTHO 2

Wavelength: 420 - 480 nm

Power: 3,000mW/cm²

Li-Ion battery - Handpieces:
anodized aluminium

Delivered with a handpiece,
a Ø5.5mm opalescent guide,
a power supply + 4 plugs,
a battery, a battery charger
and an oval light shield



MINILED — STANDARD

Wavelength: 420 - 480 nm

Power: 1,250mW/cm²

Li-Ion battery - Handpieces:
anodized aluminium

Delivered with a Mini LED
handpiece, a battery, a battery
charger, a power supply cord
with european plug,
a Ø7.5mm opalescent guide
and an oval light shield



BUY 1 x Acteon Curing Light, GET 1 x Light Guide of your choice FREE

SOFTLY Amalgamator

Softly is very stable, extremely quiet, and very easy to clean and maintain. To guarantee totally safe operation, Softly is deactivated when you open the front cover to access the capsule. Softly is made entirely of injection-molded ABS material.

It complies with all electrical and mechanical safety standards.



Contact your preferred Acteon dealer for Special Softly Pricing

Acteon Alginate Mixer

Mixes all your alginates in 5 seconds without Bubbles

Mixer features a flat digital control panel; great for applying barrier film and keeping clean. No raised buttons.

Acteon Alginate mixers provide smooth mixing for bubble free impressions. Increased tear resistance and elasticity. resulting in more accurate impressions and casts.

Less alginate wasted and more time gained.

- standardizes procedures in the dental office
- better alginate mixture will result in more accurate impressions
- easy to use, easy to clean

WARRANTY: 2 Years

Dimensions: 20x24x30vm

Alginate mixer with 4 bowls, 1 dosing bottle, 1 spatula, 1 mixer



BUY 1 x Alginate Mixer, GET 1 x Alginate Mixing Bowl FREE

All clinical applications are easy with the all-in-one air polisher

AIR N GO
EASY



- A single handpiece to cover supra-gingival prophylaxis and subgingival treatments
- Greater freedom of movement for accurate treatments
- 4 interchangeable nozzles available: SUPRA, PERIO easy, PERIO, PERIO Maintenance -
- Easier cleaning of the handpiece
- The Air-N-Go Easy handpiece is suitable for SUPRA and SUB gingival treatments.



BUY 1 x AIR N GO, GET 4 x Bottles 250g AIR-N-GO Powder FREE

(Bonus bottles contain four flavours: peppermint, lemon, raspberry & cola. o/n F10255)

EXPASYL
EXACT



An EXACT solution for atraumatic gingival retraction

Atraumatic and wider sulcular opening

Time saving 1 to 3 minutes

No risk of damage to the epithelial attachment

No pain for the patient

No bleeding/Perfect hemostasis




BUY 1 x 20 or 50 pack of Expasyl Exact GET 1 x Expasyl 5 Pack FREE

Ris|Kontrol[®] your 'must have' for safe, aseptic procedures

METAL TIPS

It is impossible to clean the inside of reusable tips adequately due to the presence of biofilm.⁽¹⁾



(1,2) Data on file

"With regard to cleaning and sterilization, disposable tips were found to be more acceptable in clinical use than metal tips."⁽¹⁾

Save time & money in your daily basis practice

- Single-use
- Avoid cleaning and sterilization steps

Ensure your clinical protection 20 times less risk means more safety and serenity on a daily basis, for safe care and greater patient satisfaction.

Protect your patient Single-use = reduced risk of direct or cross-contamination

- Minimise the spread of viruses

SAVE
75%*
WITH
RISKONTROL[®] COMPARED TO REUSABLE TIPS

DIRECT CONTAMINATION
20 TIMES LESS CROSS-CONTAMINATION RISKS⁽¹⁾
CROSS-CONTAMINATION



BUY 4 x Boxes Riskontrol tips GET 1 x Box Riskontrol FREE

The quality of the diagnosis depends on the quality of the mirror used



FRONT SURFACE

Cross-sectional view





Image reflected directly on the surface of the glass, distortion-free
NO DOUBLE IMAGE

BACK SURFACE


Cross-sectional view of an aluminium mirror



Reflecting surface under the glass
DOUBLE IMAGE

Flat and ergonomic soldering

More comfortable for the patient while retracting the cheek




Satined finition

Prevent light reflection in the eyes
Avoid eye fatigue



Optimal crimping

The frame suits the mirror
No sharp edges
No blood and dirt retention



BUY 4 x Boxes of Mirrors (Pure Reflect, Titanium or Rhodium) GET 1 x Box FREE
Buy 3 x Photographic Mirrors GET 1 x (same Coating) FREE
Buy 3 x Micro Mirrors Get 1 FREE

Complete the order form and email to acteon.anz.promo@acteongroup.com

☐ **BUY** 1 x 1 P5XS B.LED or P5 B.LED or BOOSTER
GET 1 x Handpiece **FREE**

☐ **BUY** 1 x Alginate Mixer
GET 1 x Mixing Bowl **FREE**

☐ **BUY** 1 x Piezotome Cube (Ref# F50108)
GET 1 x Extraction Tip kit II (Ref# F87546) or Universal Kit (NEW KIT) **FREE**

☐ **BUY** 1 x Implant 2 Centre (Ref# F27203)
GET 1 x Cube LED Handpiece & Cord (F12816) **FREE**

☐ **BUY** 1 x Servotome
GET 1 x Electrode **FREE**

☐ **BUY** 1 x 20 or 50 Pack Expasyl Exact
GET 1 x (5) Pack of Expasyl Exact **FREE**

☐ **BUY** 4 x Boxes Riskontol
GET 1 x Box **FREE**

☐ **BUY** 1 x Acteon Curing Light
GET 1 x Light Guide **FREE**

☐ **BUY** 4 x Boxes of Mirrors
(Pure Reflect, Titanium or Rhodium))
GET 1 x Box **FREE**

☐ **BUY** 3 x Photographic Mirrors (Same coating)
GET 1 x **FREE**

☐ **BUY** 3 x Micromirrors
GET 1 x **FREE**

☐ **BUY** 1 x Air-N-Go Easy
GET 1 x 4 Bottles 250g (4 Flavors) Air-N-Go Powder (peppermint, lemon, raspberry, cola) **FREE**

Our vision - Our mission

We believe healthcare is first and foremost about quality of life. Medical treatments are more than just procedures, they restore health, comfort, confidence, and happiness and improves life of millions of patients every year. Our vision is to pave the way of a more human-centered healthcare industry based on proactive lifetime care.

To remain at the top of their game and focus on what matters most, healthcare professionals trust integrated solutions that empower them and improve their practices on a daily basis. We provide comprehensive, intuitive workflow solutions that put technology at the service of people and help healthcare professionals restore patients health effectively, efficiently, and economically.

Please print **PROMO CODE: Con N2023 -4** legibly on your Invoice. To process your free good(s), email and attach a PDF of your proof of purchase to acteon.anz.promo@acteongroup.com

Name: _____

Office Name: _____

Account #: _____

Address: _____

Postcode: _____

Phone: _____

Email: _____

I would like to receive promotional, technical, and product information via email. Please send me an email so that I can opt-in to Acteon eNews.*

The required purchase from an authorized Acteon dealer must be reflected on one invoice dated between October 1 st to December 31 st 2023. Invoices received after January 25th 2024 will not be honoured. Please, print promo code CON N2023-4 on legibly on your invoice. To process your free good(s), the doctor will email and attach a PDF of your proof of purchase to acteon.anz.promo@acteongroup.com Multiple redemptions on a single invoice are valid; Offer cannot be combined with any other Acteon or dealer special. Allow 4 weeks for delivery of free goods. Incomplete or illegible coupons or invoices will not be honoured.. Pictures for illustration purposes only.

ACTEON
ANZ
Acteon Australia

Suite 1, 13B Narabang Way
Belrose NSW 2085

www.acteongroup.com